


Micro-XRF imaging reveals the distribution of inorganic elements inside edible crickets

Chiya NISHIMURA^[1], Yuka Yanai^[1], Shunsuke Murata^[1], Sofia GAIASCHI^[2] (¹HORIBA, Ltd., Japan, ²HORIBA FRANCE SAS, France)

1. Motivation


Elemental analysis is significant for better understanding of edible insects from nutritional and toxicological perspectives. ICP-OES and XRF have been popular elemental analysis to determine total contents of inorganic elements in insects.^[1,2] However, there is no report that discuss the distribution of the elements even though they are often eaten as they are without ground into powder.

In this poster, we introduce micro-XRF as an elemental imaging tool to visualize the distribution of inorganic elements in an insect body.



2. Sample information

a pack of edible crickets from an entomophagy shop




- Country of origin: Japan
- Feed: Almond
- Process: Roasted
- Flavoring: Nothing
- Size: approx. 30 mm

Select **2** of them to check the repeatability

3. Method

What is micro-XRF? How does it work?



Method: micro-XRF*


*It is a non-destructive elemental analysis technique based on energy dispersive X-ray fluorescence spectroscopy with a micron spot size. Thanks to the motorized XY stage, it allows distribution imaging in microscopic level.

Instrument name: XGT-9000X-ray Analytical Microscope (HORIBA, Ltd.)

Sample setting: Put them on a tape on a film*

*We fixed them on a special sample tray (hole in the center) to reduce the background signals derived from scattered primary X-rays.


Condition: 50 kV, 1000 μ A, Rh target, Partial vacuum mode



Schematic diagram of XGT-9000


4. Result

The below results show the optical images, the transmission X-ray images, and the Zn images on the crickets using a HORIBA XGT-9000.



Entire imaging covering the two crickets using 100 μ m ultra-high intensity probe (in 10 min)

Detailed imaging around a jaw using 10 μ m probe (in 60min)



Micro-XRF imaging results show **Zinc enrichment in a jaw** of both crickets. Our result was consistent with previous research reporting Zn enrichment in the jaws of grasshoppers^[3], worms^[4], and a dragonfly^[5] from the biological point of view.

*Comment: The other elements didn't have any specific enrichment in the crickets. Therefore, they are not displayed here.

5. Conclusion

Our result shows the importance of micro-XRF imaging to know if an element has any distribution in an edible insect or just exists homogeneously inside. HORIBA XGT-9000 Analytical Microscope, which has dual detectors for fluorescent X-rays and transmission X-rays, can visualize where the enrichment exists without sample pretreatment and sample destruction. It can contribute to the further spread and understanding of entomophagy.

6. Further challenges

- Collect the data of more various edible insects
- Possibility to convert the X-ray signal counts into concentration of an element

7. Reference

- [1] Ghosh et al. (2017) *J. Asia Pac. Entomol.*, 20 (2) pp.686-694.
- [2] Aydoğan et al. (2021) *Int. J. Trop. Insect Sci.* 41, pp.3049-3054.
- [3] Laiolo et al. (2021) *Evolution*, 75 (5) pp.1132-1142.
- [4] Lichtenegger et al. (2003) *Proc. Natl. Acad. Sci. USA.*, 100 (16), pp.9144-9149.
- [5] Nakamura et al. (2017) *Adv. X-Ray. Chem. Anal., Japan*, 48, pp.365-374.

